For veterinary use only

English version

SNAP* 4Dx* Plus Test
 In vitro diagnostic for the detection of *Dirofilaria immitis* antigen, antibody to *Anaplasma phagocytophylum*, antibody to *Anaplasma platys*, antibody to *Borrelia burgdorferi*, antibody to *Ehrlichia canis*, and antibody to *Ehrlichia ewingii* in canine serum, plasma, or whole blood.

Precautions and warnings

- All wastes should be properly decontaminated prior to disposal.
- Do not mix components from kits with different serial numbers.
- Do not use a SNAP device that has been activated prior to the addition of sample.
- Male-only heartworm infections typically produce antigen levels that are below the detection capability of this analyte test.
- Refer to country specific Material Safety Data Sheet for regional hazard identification.

Storage

- Store at 2–8°C until expiration date.
- SNAP devices and reagents can be stored at room temperature (18–25°C) for 90 days or until the expiration date, whichever occurs first.
- After SNAP devices and reagents are removed from 2–8°C for more than **24 hours**, the expiration date is 90 days or the printed expiration date, whichever occurs first. If the 90-day expiration date occurs prior to the printed expiration date, record the new date in the space provided on the kit.

Item	Reagents	Quantity
1	1 or 5 bottles Anti-D. <i>immitis</i> /Anaplasma spp./ <i>B. burgdorferi</i> /E. <i>canis</i> /E. <i>ewingii</i> : HRPO conjugate (preserved with gentamicin and ProCin™ 150)	7,0 mL
2	SNAP device	5, 15, 30, or 180
Reagents contained in each device:		
	Wash solution (preserved with ProCin™ 150)	0,4 mL
	Substrate solution	0,6 mL
Other components: transfer pipettes, sample tubes, and reagent rack		

Conjugate - H317/P261/P280/P302+P332/P333+P313. May cause an allergic skin reaction. Avoid breathing mist/vapors. Wear protective gloves. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.



Wash Solution - H317/H319/H402/H412/P261/P280/P305+P351+P338/P337+P313/P333+P313. May cause an allergic skin reaction. Causes serious eye irritation. Harmful to aquatic life. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapors. Wear eye protection/face protection. Wear protective gloves. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. If skin irritation or rash occurs: Get medical advice/attention.

Sample information

- Samples must be at room temperature (18–25°C) before beginning the test procedure.
- Serum, plasma, or anti-coagulated whole blood (e.g., EDTA, heparin), either fresh or stored at 2–8°C for up to one week, can be used.
- For longer storage, serum, or plasma can be frozen (20°C or colder) and then recentered before use.
- Hemolyzed or lipemic samples will not affect test results.

Test procedure

- If stored in a refrigerator, allow all components to equilibrate at room temperature (18–25°C) for 30 minutes before use. **Do not heat.**
- Using the pipette provided, dispense **3 drops of sample** into a new sample tube.
- Hold the bottle vertical, add **4 drops of conjugate** to the sample tube.
- Cap the sample tube and mix it thoroughly by **inverting it 3–5 times.**


When color **FIRST** appears in the activation circle, push the activator firmly until it is flush with the device body.

Note: Some samples may not flow to the activation circle within 60 seconds, and therefore, the circle may not turn color. In this case, press the activator after the sample has flowed across the result window.

Read test result at 8 minutes.
Note: The positive control may develop sooner, but results are not complete until 8 minutes.

Interpreting the test results
Positive result
 Any color development in the sample spots indicates the presence of *Dirofilaria immitis* antigen, *A. phagocytophylum* antibody, *A. platys* antibody, *B. burgdorferi* antibody, *E. canis* antibody, or *E. ewingii* antibody in the sample.

Negative result
 Only positive control spot develops color.

Invalid results

- Background**—If the sample is allowed to flow past the activation circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No color development**—If the positive control does not develop color, repeat the test.

IDEXX SNAPshot Dx® analyzer
 Test results can also be read using the SNAPshot Dx® analyzer. A complete description of how to enter patient data and read test results using the SNAPshot Dx analyzer can be found in the SNAPshot Dx analyzer operator's guide.

Borrelia burgdorferi vaccine cross-reactivity—The *Borrelia burgdorferi* assay detects antibodies induced as a result of natural infection with the organism and not following immunization with the following vaccines:
 Recombitek® Lyme, LymeVax®, and Galaxy® Lyme and Nobivac® Lyme.
 Recombitek is a registered trademark of Merial, Inc. LymeVax® is a registered trademark of Pfizer, Inc. Galaxy and Nobivac are registered trademarks of Merck Animal Health.

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
IDEXX **IDEXX** **IDEXX**

© 2022 IDEXX Laboratories, Inc. All rights reserved. • 06-28002-11

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