

IDEXX EliteVision™ System

Technical Specifications

IDEXX EliteVision™ System

Detector Type	Amorphous silicon on glass—no tiling
Detecting Area	9.59" x 11.69" (24.35 cm x 29.70 cm)
Plate Dimensions	H: 16.3" (415 mm) W: 11.2" (284 mm) D: 0.6" (16 mm) Wt: 4.8 lb (2.18 kg)
Imaging Area	14" x 17" (35.6 cm x 43.2 cm)
Image Resolution	139 μ m x 139 μ m
Software	CARESTREAM* Image Suite* Version 4.0
Operating/Storage Temperature	59°F–86°F (15°C–30°C)

*CARESTREAM and Image Suite are trademarks of Carestream Health, Inc.

© 2014 IDEXX Laboratories, Inc. All rights reserved. • 09-80948-00

All ®/TM marks are owned by IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.
The IDEXX Privacy Policy is available at idexx.com.