

System Requirements

- Windows® XP operating system (service pack 2) or later
- IDEXX-PACS™ Imaging Software or Capture™ Imaging Software, version 3.7.4 or later
- Internet access; preferably with a high-speed network or cable-modem connection

Recommended connection speeds

The speeds below are recommendations only. A faster Internet connection speed will improve performance.

- Image upload: 1 Mbps or faster
- Image download: 5 Mbps or faster

IDEXX recommends using an online speed test tool, such as the Speakeasy Speed Test* (<http://www.speakeasy.net/speedtest>), to determine your current upload and download speeds. Select a city near your location as the test destination.

Recommended software/hardware for image viewing

- Browsers: Chrome™, Firefox® or Internet Explorer® 8, 9 or 10 (Internet Explorer 8 may not support all features)
- Monitor with 1024x768 resolution or higher
- To view images on a tablet: Safari® browser on an iPad® tablet with iOS 5 or later, or Chrome™ browser on an Android™ tablet

*Please note: Online speed test tools are third-party web sites that are not controlled by IDEXX. IDEXX is not responsible for any content or links on a third-party site.